**READ BEFORE USING**

THE TISSUE WAS RECOVERED FROM DECEASED DONORS WHOSE LEGAL NEXT-OF-KIN HAS GIVEN PERMISSION FOR THE BONE TO BE DONATED. RECOVERY WAS PERFORMED USING ASEPTIC PROCEDURES. PROCESSING AND PACKAGING WERE PERFORMED USING ASEPTIC TECHNIQUES IN CONTROLLED CLEAN ROOM ENVIRONMENTS.

**DO NOT STERILIZE**

**Description**

Allogenix™ DBM products (Allogenix DBM Putty, Allogenix DBM Paste and Allogenix DBM Plus) contain human tissue (allograft bone) and are intended for transplantation. The allograft bone has been granulated, demineralized and provided in a lipid carrier. The lipid carrier used in the manufacture of Allogenix DBM products is extracted from soybean. The Allogenix DBM Plus products contain porous ceramic granules that are a composite of highly resorbable calcium carbonate with a slower resorbing 2 µm to 10 µm outer layer of calcium phosphate.

**Use**

Allogenix DBM products may be used for filling craniofacial defects and craniotomies that are no larger than 25 cm<sup>2</sup>. The devices are also to be used for filling bony voids or gaps in the extremities and pelvis that are not intrinsic to the bony stability of the structure. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to bone. The amount of Allogenix DBM products to be used should be based on the type of procedure and size of the graft site.

Allogenix DBM products have been processed aseptically and are ready to use. The Allogenix DBM products are provided in sterile syringes. The product is then packaged in two sterile peel pouches. The outer peel pouch can be used to aseptically transfer the sealed inner pouch to the sterile field. Follow these instructions:

1. Peel open outer pouch using proper sterile technique.
2. Pass inner pouch into sterile field.
3. Peel open inner pouch and remove product.
4. Remove cap.
5. Express from syringe.

Allogenix DBM products require no reconstitution prior to use. They do not require rehydration or any special preparation. Do not subject Allogenix DBM products or their packaging to additional disinfection or sterilization procedures. Allogenix DBM is intended for single patient use only. To prevent cross-infection of graft recipients and graft contamination, do not use the containers should be disposed of in accordance with recognized procedures for discarding medical waste materials.

**Patient Records**

The clinic or hospital is responsible for maintaining recipient records for the purpose of tracing allograft tissue post-implantation. Ensure that the following information is recorded in the patient's medical record and the hospital implant records (1-3 are required and 4-10 are suggested):

1. Description of Tissue
2. Donor Identification Number
3. Product Code
4. Expiration Date
5. Quantity Implanted
6. Antibiotics Used
7. Description of Procedure
8. Date and Time of Procedure
9. Surgeon Name
10. Any Other Pertinent Information

As a convenience, a Graft Tracing Record has been included to be completed at the time of the surgical procedure. A completed original is to be retained in the patient record and the copy sent back to Intertec Cross International, a Biomet company, as indicated on the enclosed Graft Tracing Record. If the entire tissue product was discarded, return the Graft Tracing Record and explain the reason for discard.

Once completed, return the bottom copy of the form using the self-mailer (or fax) to Intertec Cross International for our permanent records. File the top copy in the patient chart.

**Contraindications**

Allogenix DBM products are contraindicated in patients with incomplete skull growth.

Allogenix DBM products are contraindicated if active or latent infection is observed in or surrounding the implantation site.

**Warnings and Precautions**

Allogenix DBM products are aseptically processed and remain aseptic during the stated shelf-life in an unopened and undamaged package. The product must be used before the expiration date. Allogenix DBM must not be used under any of the following conditions:

- If any of the package or product elements appear to be missing, tampered with or damaged;
- If the product label or identifying bar code is severely damaged, illegible, or missing;
- If the expiration date shown on the package label has passed.

Do not subject Allogenix DBM products or their packaging to additional disinfection or sterilization procedures. Allogenix DBM is intended for single patient use only. To prevent cross-infection of graft recipients and graft contamination, do not use the contents on multiple patients. Empty or partially used containers should be disposed of in accordance with recognized procedures for discarding medical waste materials.

The production of all Allogenix DBM products is performed in environmentally controlled conditions and under rigorous quality controls. The DBM component is processed with 50 units/mL Bacitracin, 500 units/mL Polymyxin B, Isopropyl Alcohol, Hydrochloric Acid, Sodium Phosphate Buffer, Hydrogen Peroxide and Allowash® solution and may contain traces of these processing agents.

Allogenix DBM products must be appropriately placed and/or fixed according to the clinical requirements for the specific procedure to avoid potentially adverse effects.

Closed suction or drainage is recommended to prevent wound fluid accumulation when using Allogenix DBM products.

Allogenix DBM products may extrude into facial soft tissues and the effect of extrusion in cranial applications, due to lack of soft tissue, has not been investigated.

Use of Allogenix DBM may result in loss of contour.

The use of Allogenix DBM in a closed cavity may result in possible pressurization of the cavity which could produce fat embolization or embolization of device into the blood stream.

When used to fill craniofacial defects with minimal tissue coverage, care should be taken to avoid overfilling the defect site with Allogenix. Excessive pressure or tension on the overlying tissue by an overfilled defect could cause the wound to re-open. In addition, the use of Allogenix should be limited to craniofacial defects 25 cm<sup>2</sup> and smaller which are characterized by pronounced bony margins that can effectively contain the product at the site.

Each lot of DBM incorporated into Allogenix DBM is assayed for its osteoinductive potential. The assay measures the proliferation of SAOS human osteosarcoma cells in the presence of human DBM compared to positive and negative controls (osteoinductive index). Results of the assay have been correlated with results from implantation of DBM into athymic rat muscle, which demonstrated a correlation coefficient of 0.950 (p<0.0005) and accurately predicted the *in vivo* osteoinductivity in 25 donor lots.<sup>1</sup> Additionally, clinical results using DBM with >0.20 and ≤0.20 (osteoinductive index) demonstrated a significant difference in healing as evaluated by radiography, 92% and 33% healing, respectively.<sup>2</sup>

The combination of DBM, the carrier and, in some formulations, ceramic granules has not been evaluated for osteoinductivity; therefore it is unknown to what extent the formulation components may alter the osteoinductive character of the DBM. Additionally, it is unknown how osteoinductivity of the DBM component, measured via the *in-vitro* "SAOS" bioassay will correlate with human clinical performance of Allogenix DBM products.

<sup>1</sup> Addison HD, Strauss-Schoenberger J, Gillis M, Wilkins R, Jackson M, and Huskie KA. Rapid Quantitative Bioassay of Osteoinduction. J Ortho Res, 2000, 18:503-511.

<sup>2</sup> Wilkins RM. Clinical Efficacy of Demineralized Bone Matrix. *In-vitro* Human Cell Culture, Advances in Tissue Banking, 1989 3:113-124

Extensive donor blood serum testing, medical and social history screening procedures, and tissue microbiological testing have been used in the qualification of all tissue donors. Despite the extensive tissue donor selection and qualification processes used in providing this tissue graft, transmission of infectious diseases through the use of this tissue graft is still possible. Infection at the graft site may also occur. Any adverse outcomes potentially attributable to Allogenix DBM products must be reported promptly to Biomet Microfixation at (800) 874-7711 USA; (904) 741-4400

**Storage**

Allogenix DBM products should be stored in a clean, dry place at ambient temperature. Keep out of direct sunlight and do not freeze. It is the responsibility of the

tissue dispensing service and user facility/clinician to maintain any tissue in appropriate conditions prior to use.

**Production**

**Tissue Processing Procedures:**

Allogenix products are produced by a validated proprietary production process. Tissue processing is completed in accordance with procedures recommended by the American Association for Tissue Banks (AATB). The lipid carrier and ceramic granules have been irradiated with a Cobalt 60 source.

**Summary of Donor Records:**

The completed donor chart for the enclosed tissue (including but not limited to, serology results, recovery culture results, medical and social history evaluation and hemodilution calculation) that was conducted by or contract tested by and for the tissue provider) has been approved for transplantation by the tissue provider's Medical Director.

**Donor Screening and Testing:**

Prior to donation, the donor's medical/social history was screened for medical conditions or disease processes that would contraindicate the donation of tissues. The tissue provider's policies and procedures for donor screening, serologic and microbiologic testing meet current standards established by the AATB. Contraindications for tissue donation include, but are not limited to, the following: presence of infectious diseases, malignant disease, neurological degenerative diseases, disease of unknown etiology, and exposure to toxic substances. The donor's medical/social history was also screened for HIV high risk factors in accordance with current United States Public Health Services Recommendations for the Prevention of HIV Transmission through Tissue and Organ Donation.

Testing of donor blood and tissue samples began at the site of recovery and continued throughout processing. Donor blood samples taken at the time of recovery were tested and found to be negative (acceptable) for the following:

- Hepatitis B surface antigen (HBsAg)
- Hepatitis B core antibody (anti-HBc)
- Hepatitis C antibody (anti-HCV)
- Human immunodeficiency virus Type 1 and Type 2 antibodies (anti-HIV-1 and anti-HIV-2)
- Human T-lymphotropic virus Type 1 and Type 2 antibodies (anti-HTLV-1 and anti-HTLV-2)
- Syphilis (RPR or FTA)
- Nucleic Acid Testing (NAT) for HIV-1
- NAT for HCV

The individual tissues collected at recovery were subject to microbiological testing and determined to be free of specific aerobic/ anaerobic microorganisms and fungal contaminants whose presence would preclude the tissue from transplantation.

The above tests were performed by laboratories certified under the Clinical Laboratories Improvement Amendments of 1988 (CLIA).

**CAUTION**

Federal (USA) law restricts this device to sale by or on the order of a licensed physician.

1520 Tradecrest Drive  
Jacksonville, FL 32218  
Tel: (800) 874-7711 USA; (904) 741-4400  
Fax: (904) 741-4600

For additional information call  
Biomet Microfixation:  
(800) 874-7711 USA; (904) 741-4400

LD002-1670

REV. C 06/08

All trademarks are a property of Biomet, Inc. or one of its subsidiaries unless otherwise noted.

Allowash® is a registered trademark of LifeNet.

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Tüm ticari markalar, aksi belirtilmedikçe Biomet, Inc veya bağlı ortaklığımızın mülkiyetindedir.