

LactoSorb® Endobrow

Resorbable Technology



BIOMET®
MICROFIXATION
Anticipate. Innovate.™


Endobrow Lift Procedures With Proven Resorbable Technology

Anticipation and innovation. These two qualities have made Biomet Microfixation an industry leader. Founded by Walter Lorenz more than thirty years ago, Biomet Microfixation offers instrumentation, plating systems and related products for a wide range of surgical procedures.

The LactoSorb® Endobrow screw combines our push screw technology with our proven LactoSorb resorbable material to deliver an easy to use screw with predictable resorption.

Requiring only two instruments for insertion, the LactoSorb Endobrow push screws are low in profile and resorb within one year. Designed specifically for brow lift procedures, the unique screw design features an eyelet in the screw head to accept suture for securing the brow.

The LactoSorb Endobrow push screw may help make your brow lift procedures quicker and easier to perform.



Endobrow Lift Set: 915-1006

Includes

- 915-1005 LactoSorb® Endobrow Container (1)
- 915-3002 LactoSorb® large Tap Handle (2)
- 915-3007 LactoSorb® Hex Head Driver (1)
- 915-3008 1.8mm x 60mm, 5mm stop, Hand Drill (2)

Items Sold Separately

- 915-2333 2.0mm x 5.0mm Sus Push Screw (1 pack)
- 915-2336 2.0mm x 5.0mm Sus Push Screw (2 pack)
- 915-3009 1.8mm x 50mm, 5mm stop, J-Notch Drill



Pre-Operative



Post-Operative



Application Guidelines




1. Initial screw placement should be determined by the surgeon. Drill a pilot hole using the 1.8mm drill bit at the desired location for brow suspension. No tapping is required.



2. Load the Endobrow Push screw with suture passer onto the Hex Head Driver. If desired, the screw may be threaded with the suture at this step.*



3. Position the screw into the pilot hole (perpendicular to the bone) and push straight down until the screw is fully seated. No twisting is required during insertion.



4. Disengage the hex head from the screw after insertion by gently turning the Hex Head Driver clockwise.



5. Suture is easily passed through the screw head utilizing the suture passer (if screw is not already threaded with suture) and is fixed to both the periosteal tissue and the LactoSorb® Endobrow Push Screw.

*LactoSorb® Endobrow Push Screw can accept 2-0 suture without a needle and 3-0 suture with a needle

BIOMET[®]
MICROFIXATION
Anticipate. Innovate.™

What fascinates you about the body is also what drives us. That's why we're always pushing the boundaries of engineering to make products that help you keep the human form as glorious as it was intended. To learn more about our breadth of products, call 800-874-7711 or visit us online at biometmicrofixation.com. We'd love to join you in a conversation about the future.

BIOMET®
MICROFIXATION
Anticipate. Innovate.™

For more information on LactoSorb® Endobrow, please contact us at:

GLOBAL HEADQUARTERS

1520 Tradeport Drive • Jacksonville, FL 32218-2480
Tel (904) 741-4400 • Toll-Free (800) 874-7711 • Fax (904) 741-4500 • Order Fax (904) 741-3059
www.biometmicrofixation.com

EUROPE

Toermalijnring 600 • 3316 LC Dordrecht • The Netherlands
Tel +31 78 629 29 10 • Fax 31 78 629 29 12
e-mail: europa@biomet.com

As the manufacturer of this device, Biomet Microfixation does not practice medicine and does not recommend this product or surgical technique for use on a specific patient. The surgeon who performs any implant procedure must determine the appropriate device and surgical procedure for each individual patient. Devices shown in this brochure may not be cleared or licensed for use or sale in your individual country. Please contact your local distributor for information regarding availability of this product. Information contained in this brochure is intended for surgeon or distributor information only and is not intended for patient distribution. All surgeries carry risks. For additional information, please see package insert, or visit our web site at www.biometmicrofixation.com or call 1-800-874-7711. LactoSorb® is a registered trademark of Biomet, Inc.