BIOMET® MICROFIXATION

Biomet Microfixation
1520 Tradeport Drive
Jacksonville FL 32218 USA

LD02-4150

Rev. A 6/08

Allogenix™ LT Demineralized Bone Matrix (1cc)

ATTENTION OPERATING SURGEON

DESCRIPTION

Allogenix™ LT Demineralized Bone Matrix (Allogenix™ LT) is processed human bone that has been demineralized and combined with human collagen-derived carrier from the same donor. The final demineralized bone matrix is in a freeze-dried state. Allogenix™ LT is supplied in single-use packages for single-patient use.

Allogenix™ LT contains donated human tissue procured from human cadaveric donors. The tissue has been determined eligible for transplantation by a qualified tissue bank medical director after review of medical and social history, hospital records, infectious disease screening, autopsy report (if performed), and physical exam. Donors are tested and found negative (acceptable) for anti-HIV 1/2, HBsAg, anti-HBc, anti-HCV, anti-HTLV I / II, syphilis, HCV NAT, and HIV NAT U.S. Food and Drug Administration (FDA) licensed test kits are used when available for a specific test. Communicable disease testing has been performed by a laboratory certified under CLIA or equivalent requirements.

Before demineralization, the tissue has been processed with Allowash®, a patented bone and soft tissue cleaning technology under license from LifeNet.

Allogenix™ LT is processed and prepared via a proprietary process at Interpore Cross International, Irvine, California.

INDICATIONS FOR USE

Allogenix™ LT can be used to fill bony voids or gaps that have been surgically created, or for filling osseous defects in non-weight bearing applications.

Allogenix™ LT may be used with orthopedic, spinal, reconstructive, craniofacial, maxillofacial and periodontal bone grafting procedures. It can be used alone or in combination with autologous bone, or other forms of allogenic bone in grafting procedures of non-weight bearing value. It can be used or hydrated with autologous blood, bone marrow aspirate or, autologous blood derived products such as platelet rich plasma and platelet poor plasma. It may also be hydrated with saline or antibiotic solution.

CONTRAINDICATIONS

Infection at the surgical site and/or distant foci of infections that may spread to the surgical site is a contraindication for Allogenix™ LT.

RELATIVE CONTRAINDICATIONS

1. Uncooperative patient, or patient with neurologic disorders who is incapable of following directions, including weight control and activity levels.
2. Pregnancy.
3. Disorders or diseases which may impair bone formation.

WARNINGS AND PRECAUTIONS

Patient selection factors to be considered should include: 1. the ability and willingness of the patient to follow instructions; 2. control of weight and activity levels; and /or 3. a good nutritional state.

1. Allogenix™ LT contains donated human tissue.
2. This tissue has been processed with Bacitracin and/or Polymyxin B, HCl, alcohol, and sodium phosphate. Traces may remain.
3. Although this tissue has been tested and screened for selected human pathogens, processed under aseptic conditions, and gamma irradiated with a Cobalt 60 source at 15-25 kGy, human derived tissue may still transmit infectious agents.
4. Do not use Allogenix™ LT if package integrity has been compromised.
5. This tissue is intended for use in one patient on a single occasion only.
6. Once user breaks the container seal, the tissue must be transplanted or discarded.
7. This tissue may not be sterilized or re-sterilized.
8. Do not use for treatment of bone with compromised stability or load bearing value, or within articulating joints.

9. The surgeon is to be thoroughly familiar with Allogenix™ LT material and the surgical procedure prior to use of this tissue.
10. The patient is to be made aware of general risks associated with treatment and the possible adverse effects.
11. The product must not be used if the expiration date shown on the package label has passed.

POSSIBLE ADVERSE EFFECTS

1. Complications associated with surgery such as hematoma, infection, migration, and other complications that may require additional surgery.
2. Incomplete or lack of bony ingrowth at the treatment site that may require additional surgery.
3. Immune rejection of the introduced tissue that may require additional surgery.
4. The transmission of known pathogens including Human Immunodeficiency Virus 1/2, Hepatitis B and C, Human T-Lymphotropic Virus I and II, Syphilis, bacteria, and fungi.

STERILITY

This tissue has been processed under aseptic conditions. Allogenix™ LT has been irradiated in its final container with a Cobalt 60 source at 15-25 kGy.

INSTRUCTIONS FOR USE

1. Twist the elbow to seat perpendicular to the tabs of the 1cc Graft Preparation System containing Allogenix™ LT. Ensure that the elbow attachment is well seated in the 1cc syringe with Allogenix™ LT. If not, push on elbow to seat (Figure 1).
2. Ensure that the yellow cap on the Graft Preparation System containing Allogenix™ LT is tightly seated.
3. Attach the hydrating fluid of choice by twisting the syringe onto the Luer Lock fitting. Ensure a minimum ratio of 0.6ml fluid to 1cc DBM prior to attaching syringe to Graft Preparation System. **NOTE** – in all subsequent steps, the syringe with hydrating fluid must be pointing down as illustrated in Figure 2.
4. Pull on the syringe plunger containing the hydrating fluid until fully extended. This will create a vacuum necessary for complete hydration (Figure 3).
5. Release plunger to allow the hydrating fluid to infiltrate the Graft Preparation System. Wait 2 to 3 minutes before proceeding (Figure 4).
6. Use hydrating syringe to lever out the elbow. Attach plunger provided, remove yellow cap and expel hydrated Allogenix™ LT (Figure 5).

STORAGE AND SHELF LIFE

Ambient temperature is recommended. No refrigeration is necessary. See package label for expiration date. It is the responsibility of the Tissue Dispensing Service and/or end-user to maintain Allogenix™ LT in the appropriate storage conditions prior to transplant.

TRACKING AND TRACEABILITY

Please complete the enclosed Graft Tracking Record and return it to Interpore Cross International, a Biomet company, following the directions provided on the Graft Tracking Record. Federal regulations (21 CFR.290 (b)) require proper tracking of human tissue. It is the responsibility of the end-user to provide this tracking information, which enables Biomet Microfixation to maintain records for the purpose of human tissue post-transplant or any other final disposition (e.g. tissue not used and discarded). Adverse outcomes potentially attributable to the tissue must be promptly reported to Biomet Microfixation, Inc. Use the peel-off sticker from the label in the patient records.

Caution: Federal law (USA) restricts this tissue to sale, distribution, or use by or on the order of a licensed physician. This tissue is intended for use by qualified health care specialists such as physicians, dentists, or podiatrists.

Biomet Microfixation and Interpore Cross International make no claims concerning the biological or biomechanical properties of the provided tissue. All tissue is recovered, processed, stored, and distributed for use in accordance with the standards of the American Association of Tissue Banks. Biomet Microfixation and Interpore Cross International disclaim all liability and responsibility for any misuse of tissue provided for clinical application.

Comments regarding this tissue can be directed to Attn: Regulatory Dept., Biomet Microfixation 1520 Tradeport Dr. Jacksonville FL 32218 USA Phone: (800) 874-7711 USA; (904) 741-4400 FAX: (904) 741-4500 Patent No. 6,576,249 and other pending patents

Biomet® and all other trademarks herein are the property of Biomet, Inc. or its subsidiaries.

Figure 1


Figure 2


Figure 3


Figure 4


Figure 5

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Biomet® ve burada bahsi geçen diğer tüm ticari markalar Biomet, Inc. ya da yan kuruluşlarının mülkiyetindedir.