



The Lorenz Office System is a fully functional titanium system in a convenient container that can be tailored to fit individual surgeon's needs. This system proves that convenience and functionality do not have to be sacrificed for size. Features of the Lorenz Office System include:

- ▶ **Convenient Size and Design:** Ideal size for any office autoclave, transparent lid allows full view of inventory.
- ▶ **Fully Customizable:** System components are customized to individual office needs.
- ▶ **Flexible Storage Capabilities:** Includes a large generic storage bay for plates, mesh, instruments, grafts, distraction devices, and dental implants.
- ▶ **Removable Screw Dial:** The portable screw dial allows easy access to screws. Interchangeable size indicators and easy reload makes inventory management simple.
- ▶ **IMF Screws:** An innovative addition to the office set, these screws allow the surgeon to use IMF in the office setting.
- ▶ **Instrumentation:** Options include: Screwdriver handle and blades, membrane punch, graft fixation forceps, wire/mesh cutters, membrane forceps and scissors, drill bits, screws and other implants.
- ▶ **Affordability:** One low price includes instrumentation and implants.



01-6851

1.5mm IMF Set:

Container, dial, lid, wire cutters, forceps, screwdriver handle

01-6861

1.5/2.0mm Set

Container, dial, lid, wire cutters, forceps, screwdriver handle

Tray Components

Item #	Description
01-6850	1.5/IMF Tray, dial, and lid
01-6860	1.5/2.0 Tray, dial, and lid
01-6852	1.5/IMF Dial
01-6862	1.5/2.0 Dial
01-6858	Plate Insert Tray

Instrumentation

Item #	Description
51-1362	Wire Cutters, TC Serrated
01-6958	Membrane Punch Forceps
01-7390	Green Screwdriver Handle
01-7164	Silver Screwdriver Handle
51-0240	Iris Membrane Scissors
50-4080	Small Contra Angle Handle
51-9761	Adson Membrane Forceps

Screws

1.5 HT, X-Drive

Item #	Description
91-1502	X 2.5mm
91-1503	X 3.5mm
91-1504	X 4mm
91-1505	X 5mm
91-1506	X 6mm
91-1507	X 7mm
91-1508	X 8mm
91-1509	X 9mm
91-1511	X 11mm
91-1513	X 13mm
91-1515	X 15mm

Emergency 1.8 X-Drive

91-1803	X 3.5mm
91-1805	X 5mm
91-1807	X 7mm

1.5 HT Self Drilling

91-6103	X 3.5mm
91-6104	X 4mm
91-6105	X 5mm
91-6106	X 6mm
91-6107	X 7mm

1.5 Center-Drive

Item #	Description
99-7202	X 2.5mm
99-7203	X 3.5mm
99-7204	X 4mm
99-7205	X 5mm
99-7206	X 6mm
99-7207	X 7mm
99-7208	X 8mm
99-7209	X 9mm
99-7211	X 11mm
99-7213	X 13mm
99-7215	X 15mm

Emergency 1.8 Center-Drive

99-7223	X 3.5mm
99-7225	X 5mm
99-7227	X 7mm

1.5 HT Center-Drive Self Drilling

99-7283	X 3.5mm
99-7284	X 4mm
99-7285	X 5mm
99-7296	X 6mm

2.0 HT, X-Drive

Item #	Description
91-2005	X 5 mm
91-2007	X 7 mm
91-2009	X 9 mm
91-2011	X 11 mm
91-2013	X 13 mm
91-2015	X 15 mm
91-2017	X 17 mm

Emergency 2.3 HT, X-Drive

91-2305	X 5 mm
91-2307	X 7 mm
91-2309	X 9 mm
91-2311	X 11 mm

2.0 HT Self-Drilling

91-6204	X 4 mm
91-6205	X 5 mm
91-6207	X 7 mm

2.0 Center-Drive

Item #	Description
99-9474	X 4mm
99-9475	X 5mm
99-9477	X 7mm
99-9479	X 9mm
99-9481	X 11mm
99-9483	X 13mm
99-9485	X 15mm
99-9487	X 17mm
99-9489	X 19mm

Emergency 2.3 Center-Drive

99-9495	X 5mm
99-9497	X 7mm
99-9499	X 9mm

2.0 IMF Screws

91-5205	X 5mm
91-5207	X 7mm
91-5209	X 9mm
91-5211	X 11mm
91-5213	X 13mm
91-5215	X 15mm

Drills, Blades

1.1 Drills

Item #	Description
01-7140	1.1 X 50mm
01-7142	1.1 X 50mm, 5mm Stp
01-7144	1.1 X 50mm, 7mm Stp
01-7148	1.1 X 50mm, 15mm Stp

1.5 Drills

01-9193	1.5 X 50mm, 7mm Stp
01-9195	1.5 X 50mm, 10mm Stp
01-9196	1.5 X 50mm, 50mm No Stp

Contra Angle

Item #	Description
01-6950	1.1 X 8mm Contra Drill
01-6951	1.1 X 16mm Contra Drill
01-7162	1.5 X-Lock Blade Contra
01-7693	2.0 X-Lock Blade Contra
01-9180	1.5 X 7mm Contra Drill
01-9181	1.5 X 15mm Contra Drill
01-9182	1.5 X 18mm Contra Drill

Blades

Item #	Description
15-1196	1.5 HT Blade Std
15-1194	1.5 HT Blade Short
20-1194	2.0 X-Lock Blade Short
20-1196	2.0 X-Lock Std
20-1199	IMF Sleeve

